


Shutdown Procedure for all Single Quads and Triple Quads



Quick Shutdown of TQ & SQ running MassHunter or CS OpenLAB

Quick Shutdown (active time 3 min, total time 25 min)

- ✓ Click off
- ✓ Close ballast valve (if you have an Edwards rough pump)
- ✓ Click Vent
- ✓ Wait 20min
- ✓ Verify turbo speeds
- ✓ Turn off MS and LC modules
- ✓ Close gas supply lines to MS
- ✓ Unplug the MS power cord

Shutdown Procedure

The screenshot shows the Agilent MassHunter Workstation Data Acquisition software interface. The interface is divided into several windows, each with a label and a blue arrow pointing to it:

- Actuals window:** A table showing instrument parameters and their values. The table has two columns: "Parameter" and "Value".
- Instrument Status window:** A panel showing the status of various instrument components: HIP Sampler (Idle), Quat. Pump (Standby), Column Comp. (Not Ready), and QQQ (Standby).
- Chromatogram Plot window:** A plot showing Total Ion Chromatogram (TIC) with a peak at approximately 124.01 minutes.
- Spectrum Pane window:** A mass spectrum plot showing relative intensity versus m/z (mass-to-charge ratio).
- Method Editor window:** A panel for editing acquisition methods, including sections for "Acquisition" and "Scan parameters".

Parameter	Value
QQQ: Ready Type	instrument_state_temperature/flow
QQQ: Ready State	false
QQQ: Run Time	0.00 min
QQQ: Run State	standby
QQQ: Rough Vac	2.42E-0 Torr
QQQ: Tube 1 Speed	100.0%
QQQ: Tube 2 Speed	100.0%
QQQ: High Vac	1.04E-5 Torr
QQQ: Instrument Stage	standby
QQQ: Film State	no film

Segment Name	Start Mass	End Mass	Scan Time	Frequency	Cell Accelerator Voltage	Polarity
1	100	1000	500	300	5	Positive

#	Start Time	Scan Type	Div Valve	Delta EHV (eV)	Delta EHV (V)	Stored
1		MS2 Scan	To MS	500	500	<input checked="" type="checkbox"/>

Adding Vacuum parameters to Actuals panel

Parameter	Value
QQQ: Not Ready Text Long	
QQQ: Run State	prerun
QQQ: High Vac	1.27E-5 Torr
QQQ: Gas Temp	350 °C
QQQ: Rough Vac	1.66E+0 Torr
QQQ: Turbo 1 Speed	100.0 %
QQQ: Turbo 2 Speed	Setup... 0.0 %
QQQ: Firmware Version	A.00.07.03
QQQ: Sheath Gas Flow	12.0 l/min
QQQ: Sheath Gas Temp	400 °C
Binary Pump: Ripple	0.22 %
QQQ: Instrument State	background_acquisition

Actuals Selection Dialog

Available Parameters

- HIP Sampler
- Binary Pump
- Column Comp.
- DAD
- QQQ
 - APPI Lamp
 - Calibrant
 - Capillary
 - Capillary Current
 - Chamber Current
 - Collision Gas
 - Corona Voltage
 - Diverter Valve
 - Error State
 - Firmware Version
 - Gas Flow
 - Gas Temp
 - High Vac
 - Instrument State
 - Ion Mode
 - Ion Polarity
 - Ion Source
 - Locked Mode
 - MS 1 Heater
 - MS 2 Heater
 - Nebulizer
 - Not Ready Text Long
 - Not Ready Text Short
 - Purge 1 Current

Parameters to display

- QQQ: Not Ready Text Long
- QQQ: Run State
- QQQ: High Vac
- QQQ: Gas Temp
- QQQ: Rough Vac
- QQQ: Turbo 1 Speed
- QQQ: Turbo 2 Speed
- QQQ: Firmware Version
- QQQ: Sheath Gas Flow
- QQQ: Sheath Gas Temp
- Binary Pump: Ripple


Buttons: Clear All, Up, Down, Move To Top

Color Formatting: Background Color, Text Color, Reset

Conditional Formatting for Selected Parameter Values: Lower Limit, Upper Limit

Buttons: OK, Cancel

Close Ballast valve (if you have an Edwards rough pump)



MS40 pump (internally ballasted)



Click Off, Right click on QQQ

The screenshot shows the 'Instrument Status' window with the following components and their status:

- HIP Sampler:** Idle, EMF ✓, 5.00 µL
- Quat. Pump:** Standby, EMF ✓, 100.0 (A), 0.0 (B), 0.000 mL/min, 1.23 bar, 0.0 (C), 0.0 (D)
- Column Comp.:** Idle, EMF ✓, 21.60 °C, 22.85 °C
- DAD:** Idle, EMF ✓
- QQQ:** Idle, EMF ✓, AJS ESI

At the bottom right, the status is 'Instrument Standby' with an 'Off' button highlighted by a red box and a red arrow labeled '1'. A red arrow labeled '2' points to the QQQ component with the text '2. Right-Click'.

Click Vent

The screenshot displays the 'Instrument Status' window with five main components: Multisampler, Binary Pump, Column Comp., VWD, and QQQ. The Multisampler is in a 'Not Ready' state, while the others are 'Idle'. A context menu is open over the QQQ component, listing various control options. The 'Vent' option is highlighted with a red rectangular box.

Component	Status	EMF	Parameters
Multisampler	Not Ready	EMF ✓	0.00 µL
Binary Pump	Idle	EMF ⚠	80.00, 20.00, 0.000 mL/min, 0.00 bar
Column Comp.	Idle	EMF ✓	30.36°C, 31.30°C
VWD	Idle	EMF ✓	250 nm
QQQ	Idle	On/Off	AJS ESI, On, StandBy, Calibrant, Collision Gas, LC, Vent, Pump Down, APPI UV Lamp

Wait 20 minutes– Verify venting



Parameter	Value
QQQ: Not Ready Text Long	
QQQ: Run State	prerun
QQQ: High Vac	1.27E-5 Torr
QQQ: Gas Temp	350 °C
QQQ: Rough Vac	1.66E+0 Torr
QQQ: Turbo 1 Speed	100.0 %
QQQ: Turbo 2 Speed	Setup... 0.0 %
QQQ: Firmware Version	A.00.07.03
QQQ: Sheath Gas Flow	12.0 l/min
QQQ: Sheath Gas Temp	400 °C
Binary Pump: Ripple	0.22 %
QQQ: Instrument State	background_acquisition

Rough Pump will turn off when Turbo 1 Speed is < 10%.

Turn off power after venting



1. Turn off Power Switch **first**



2. Turn off Pump Expander box*


Installed on: **6420, 6460C, 6470, 6490, 6495**

Quick Shutdown of SQ ChemStation Edition


Quick Shutdown (active time 3 min, total time 25 min)

- ✓ Click off
- ✓ Close ballast valve (if you have an Edwards rough pump)
- ✓ Click Diagnosis Tab
- ✓ Click Maintenance Tab
- ✓ Click Vent MSD (follow steps on the screen)
- ✓ Turn off MS and LC modules
- ✓ Close gas supply lines to MS
- ✓ Unplug the MS power cord

Vent- Diagnosis Menu- ChemStation



Vent- Maintenance tab



Vent

Vent

