

Instructions for Enabling Compound Math on B.03.0x-B.07.00

Note: A new concept of Compound Math was introduced with MassHunter Quantitative Analysis B.07.01 making the calculations straightforward. Several new videos are available to introduce and explain the parameter setup.

Navigate to the following directory \Program Files\Agilent\MassHunter\Workstation\Quant\bin

- 1) Make a backup copy of the following file:

Component.config.xml

- 2) Edit the original component.config.xml with Notepad.
- 3) Uncomment the compound math portion. Before and after images shown below.

```

component.config.xml.bak - Notepad
File Edit Format View Help
<?xml version="1.0" encoding="utf-8"?>
<components>
  <defaultnamespace>Agilent.MassSpectrometry.DataAnalysis.Quantitative</defaultnamespace>
  <references>
    <reference framework="true" assembly="System.Data.dll"/>
    <reference framework="true" assembly="System.Xml.dll"/>
    <reference framework="true" assembly="System.Web.Services.dll"/>
    <reference framework="false" assembly="Quantitation.dll"/>
    <reference framework="false" assembly="Quantutils.dll"/>
    <reference framework="false" assembly="ScriptableInterfaces.dll"/>
  </references>
  <!--
  <interface name="IScriptableSample">
    <implementation assembly="CompoundMath_Script.dll">
      <script name="CompoundMath_Script.txt">
        <references>
          <reference framework="false" assembly="BatchAnalysis.dll"/>
        </references>
      </script>
    </implementation>
  </interface>
  -->
  <interface name="ICustomExpressions">

```

Agilent Technologies, Inc. believes this document to be accurate and up to date, but will not assume responsibility for its use; individuals must determine its fitness for their purposes. Microsoft, MS-DOS, .NET, Windows Vista, and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Should look as follows:

```

component.config.xml - Notepad
File Edit Format View Help
<?xml version="1.0" encoding="utf-8"?>
<components>
  <defaultnamespace>Agilent.MassSpectrometry.DataAnalysis.Quantitative</defaultnamespace>
  <references>
    <reference framework="true" assembly="System.Data.dll"/>
    <reference framework="true" assembly="System.Xml.dll"/>
    <reference framework="true" assembly="System.Web.Services.dll"/>
    <reference framework="false" assembly="Quantitation.dll"/>
    <reference framework="false" assembly="Quantutils.dll"/>
    <reference framework="false" assembly="ScriptableInterfaces.dll"/>
  </references>
  <interface name="IScriptableSample">
    <implementation assembly="CompoundMath_script.dll">
      <script name="CompoundMath_Script.txt">
        <references>
          <reference framework="false" assembly="BatchAnalysis.dll"/>
        </references>
      </script>
    </implementation>
  </interface>
  <!--
  <interface name="ICustomExpressions">

```

4) Save the file.


Next portion – add a new compound with math as part of the name.

temazepam	4	301.0 -> 255.0	MRM	Target	301.0	255.0	11.594
triazolam	4	343.0 -> 308.0	MRM	Target	343.0	308.0	10.521
diazepam	5	285.0 -> 222.0	MRM	Target	285.0	222.0	13.616
diazepam+triazolam	5	285.0 -> 222.0	MRM	Target	285.0	222.0	13.616
diazepam/triazolam	5	285.0 -> 222.0	MRM	Target	285.0	222.0	13.616

5) In the example above, input diazepam/triazolam. The transition is unimportant as the “math compound” is calculated from the two individual results, but it’s reasonable to put in an actual working value.

Agilent Technologies, Inc. believes this document to be accurate and up to date, but will not assume responsibility for its use; individuals must determine its fitness for their purposes. Microsoft, MS-DOS, .NET, Windows Vista, and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.


- 6) The “math compound” does require a calibration table. Copy a calibration table from any existing target compound to the “math compound”.



- 7) Exit and apply the method to the batch.
8) Analyze the batch.

Agilent Technologies, Inc. believes this document to be accurate and up to date, but will not assume responsibility for its use; individuals must determine its fitness for their purposes. Microsoft, MS-DOS, .NET, Windows Vista, and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

- 9) Add a column to see the custom calculation results. Note that the custom calculation is displayed in the “Compound Method” and not the “Compound Results”.



- 10) View the results. Notice that for the compound “diazepam/triazolam” the Custom Calc. shows a number close to 1, as would be expected for the calibrators.

5: diazepam/triazolam		ISTD:
diazepam/triazolam Method		
Date-Time	Exp. Conc.	Custom Calc.
107 1:24 PM		0.6217
107 1:50 PM	10.0000	0.9215
107 2:15 PM	25.0000	1.0475
107 2:41 PM	50.0000	1.0473
107 3:07 PM	125.0000	1.0168
107 3:33 PM	250.0000	0.9893
107 3:59 PM	500.0000	0.9966
107 4:25 PM		0.9617
107 4:50 PM		4.1469

Agilent Technologies, Inc. believes this document to be accurate and up to date, but will not assume responsibility for its use; individuals must determine its fitness for their purposes. Microsoft, MS-DOS, .NET, Windows Vista, and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.