



Agilent AssayMAP Bravo Cartridges for Automated Protein Sample Preparation Workflows

Selection Guide

Introduction

The Agilent AssayMAP platform is a powerful yet easy-to-use automation solution specifically designed for biomolecule sample preparation. AssayMAP microchromatography cartridges combined with optimized protocols enable non-automation experts to use automation to increase precision and scalability while at the same time decreasing the labor and human error so you have more confidence in your results and more time to analyze your biomolecule. The AssayMAP platform supports a broad range of protein quantification and characterization workflows including affinity purification, enzymatic digestion, protein and peptide cleanup, peptide mapping, peptide fractionation, N-Glycan analysis, and phosphopeptide enrichment.



Agilent Technologies

Cartridge racks

Product	Description	Application	Catalog number
C18 cartridges	1 rack (96 cartridges)	Desalting mixtures of peptides and small proteins prior to analysis by mass spectrometry.	5190-6532
Fe(III)-NTA cartridges	1 rack (96 cartridges)	Phosphopeptide enrichment prior to analysis by mass spectrometry.	G5496-60085
Protein A (PA-W) cartridges	1 rack (96 cartridges)	Purification of monoclonal and polyclonal antibodies and Fc fusion proteins with affinity for Protein A.	G5496-60000
Protein G (PG-W) cartridges	1 rack (96 cartridges)	Purification of monoclonal and polyclonal antibodies and Fc fusion proteins with affinity for Protein G.	G5496-60008
Resin-Free cartridge	1 rack (96 cartridges)	Placeholder when working with sample numbers that are not a multiple of eight.	G5496-60009
Reversed Phase (RP-S) cartridges	1 rack (96 cartridges)	Desalting or fractionating mixtures of peptides and small proteins prior to analysis by mass spectrometry.	G5496-60033
Reverse Phase (RP-W) cartridges	1 rack (96 cartridges)	Desalting of denatured proteins (antibodies) prior to trypsin digestion.	G5496-60086
Streptavidin (SA-W) cartridges	1 rack (96 cartridges)	Generation of affinity cartridges by immobilizing biotinylated molecules. Affinity cartridges can subsequently be used to purify target molecules from complex matrices.	G5496-60010
Strong Cation Exchange (SCX) cartridges	1 rack (96 cartridges)	Fractionating mixtures of peptides and small proteins prior to analysis by mass spectrometry.	5190-6536
TiO2 cartridges	1 rack (96 cartridges)	Phosphopeptide enrichment prior to analysis by mass spectrometry.	G5496-60016

Starter kits (including cartridges and required labware)

Product	Description	Application	Catalog number
AssayMAP Digestion and C18 Cleanup Starter Kit	1 rack of C18 cartridges plus labware for digestion and cleanup protocols	Protein digestion and subsequent peptide cleanup prior to analysis by mass spectrometry.	G5496-60013
AssayMAP Digestion and RP-S Cleanup Starter Kit	1 rack of RP-S cartridges plus labware for digestion and cleanup protocols	Protein digestion and subsequent peptide cleanup prior to analysis by mass spectrometry.	G5496-60034
AssayMAP SCX Fractionation Starter Kit	1 rack of SCX cartridges plus labware for digestion and cleanup protocols	Protein digestion and subsequent peptide fractionation prior to analysis by mass spectrometry.	G5496-60014

Supplemental items

Product	Description	Application	Catalog number
Syringe Test Kit	96 resin free cartridges (1 rack plus enough syringe test solution to perform 20 tests)	Routine verification that all 96 the AssayMAP syringes plus testing solution are functioning properly.	G5496-60050
Syringe Replacement Kit	10 replacement syringes plus the tools required to replace faulty syringes.	Allows customers to replace the syringes on site to avoid costly down time.	G5409-68002
AssayMAP Syringes	10 AssayMAP syringes	AssayMAP syringes for use with the syringe replacement kit.	G5496-60007
Bravo Plate Risers	1 Plate Riser	Accessory required for normalization protocol.	G5498B#061
96 well Greiner U,V-bottom plate thermal insert	Thermal insert for the Peltier heater/cooler to improve heat transfer to the plate wells.	For In-Solution Digestion Single plate and On-Cartridge Reaction applications.	G5498B#126
96 well Abgene deep well plate thermal insert	Thermal insert for the Peltier heater/cooler to improve heat transfer to the plate wells.	For In-Solution Digestion Single plate and On-Cartridge Reaction applications.	G5498B#127
Bravo closed corner plate pad	Plate pad designed to minimize the risk of plate misalignment.	Alternative for all standard Bravo plate pads.	G5498B#125
PCR Plate Thermal Insert	Thermal insert for the Peltier heater/cooler to improve heat transfer to the plate wells.	For In-Solution Digestion Single plate and On-Cartridge Reaction applications.	G5498B#013

Additional Cartridge Details

Cartridge		
C18		
Part number	5190-6532	
Bead type	Silica (20 µm diameter)	
Pore size	150 Å	
Surface chemistry	C18 (reversed-phase)	
Dynamic binding capacity	400 µg human insulin (≤10% breakthrough when loading up to 400 µg human insulin per cartridge) Approximately 150–200 µg of digested BSA (loading more than 150 µg of digested BSA results in loss of very hydrophilic peptides)	
Saturation binding capacity	Approximately 450 µg human insulin	
Storage temperature	Room temperature	
pH stability	2–7 for routine use	

Cartridge		
Fe(III)-NTA		
Part number	G5496-60085	
Bead type	Synthetic Polymer (50 µm diameter)	
Pore size	Large enough for large biomolecules	
Surface chemistry	Fe(III) chelated to nitrilotriacetic acid	
Dynamic binding capacity	20 µg adenosine monophosphate (≥ 90% recovery when 20 µg of adenosine monophosphate is loaded per cartridge) Up to approximately 800 µg of digested yeast extract can be loaded per cartridge with minimal phosphopeptide breakthrough	
Saturation binding capacity	Approximately 25 µg adenosine monophosphate	
Storage temperature	Room temperature	
pH stability	2–11 (cartridge will turn more yellow/gold with increasing pH)	
PA-W		
Part number	G5496-60000	
Bead type	Synthetic polymer (53–61 µm diameter)	
Pore size	Large enough for large biomolecules such as antibodies	
Surface chemistry	Recombinant Protein A	
Dynamic binding capacity	100 µg human IgG1 (≤10% breakthrough when loading up to 100 µg human IgG1 per cartridge)	
Saturation binding capacity	Approximately 125 µg hIgG1	
Storage temperature	Room temperature	
pH stability	6–8 for routine use 3–11 for short exposure	
PG-W		
Part number	G5496-60008	
Bead type	Synthetic polymer (45–53 µm diameter)	
Pore size	Large enough for large biomolecules such as antibodies	
Surface chemistry	Recombinant Protein G	
Dynamic binding capacity	100 µg human IgG (≤10% breakthrough when loading up to 100 µg human IgG per cartridge)	
Saturation binding capacity	Approximately 150 µg IgG	
Storage temperature	Room temperature	
pH stability	6–8 for routine use 3–11 for short exposure	
Resin Free		
Part number	G5496-60009	
Bead type	NA	
Pore size	NA	
Surface chemistry	NA	
Dynamic binding capacity	NA	
Saturation binding capacity	NA	
Storage temperature	Room temperature	
pH stability	NA	
RP-S		
Part number	G5496-60033	
Bead type	Synthetic polymer (15–20 µm diameter)	
Pore size	100 Å	
Surface chemistry	Underivatized polystyrene-divinylbenzene (reversed-phase)	
Dynamic binding capacity	400 µg human insulin (≤10% breakthrough when loading up to 400 µg human insulin per cartridge) Approximately 150–200 µg of digested BSA (loading more than 150 µg of digested BSA results in loss of very hydrophilic peptides)	
Saturation binding capacity	Approximately 500 µg human insulin	
Storage temperature	Room temperature	
pH stability	1–14	

Cartridge		
RP-W		
Part number	G5496-60086	
Bead type	Synthetic polymer (30 µm diameter)	
Pore size	Large enough for large biomolecules such as antibodies	
Surface chemistry	Underivatized polystyrene-divinylbenzene (reversed-phase)	
Dynamic binding capacity	75 µg denatured human IgG1 (≤10% breakthrough when loading up to 75 µg denatured human IgG1 per cartridge)	
Saturation binding capacity	Approximately 90 µg denatured human IgG1	
Storage temperature	Room temperature	
pH stability	1–14	
SA-W		
Part number	G5496-60010	
Bead type	Synthetic polymer (45–53 µm diameter)	
Pore size	Large enough for large biomolecules such as antibodies	
Surface chemistry	Recombinant Streptavidin	
Dynamic binding capacity	50 µg biotinylated BSA (≤10% breakthrough when loading up to 50 µg biotinylated BSA per cartridge) Approximately 90 µg biotinylated IgG	
Saturation binding capacity	Approximately 60 µg biotinylated BSA Approximately 100 µg biotinylated IgG	
Storage temperature	4–10°C	
pH stability	6–8 for routine use 3–11 for short exposure	
SCX		
Part number	5190-6536	
Bead type	Silica (20 µm diameter)	
Pore size	300 Å	
Surface chemistry	Polysulfoethyl A™ (Strong Cation Exchange)	
Dynamic binding capacity	800 µg human insulin (≤10% breakthrough when loading up to 800 µg human insulin per cartridge) Approximately 400-500 µg of digested BSA	
Saturation binding capacity	Approximately 1,000 µg human insulin	
Storage temperature	Room temperature	
pH stability	2–7 for routine use	
TiO₂		
Part number	G5496-60016	
Bead type	Ceramic titanium dioxide (20 µm diameter)	
Pore size	100 Å	
Surface chemistry	Titanium dioxide	
Dynamic binding capacity	80 µg phenyl phosphate (≥ 90% recovery when 80 µg of phenyl phosphate is loaded per cartridge) Up to approximately 1,000 µg of digested α-casein can be loaded per cartridge with minimal phosphopeptide breakthrough	
Saturation binding capacity	Approximately 115 µg phenyl phosphate	
Storage temperature	Room temperature	
pH stability	1–14	

For Research Use Only. Not for use in diagnostic procedures. This information is subject to change without notice.

©Agilent Technologies, Inc., 2015, 2016, 2019
Published in the USA, January 22, 2019
5991-4863EN



Agilent Technologies