

G1312B**Pressure Test****Aborted by system**

DEADS10554
A.07.01[001]
1260 Bin Pump

The pressure test is a quick, built-in test designed to demonstrate the pressure-tightness of the system. The test should be used when problems with small leaks are suspected, or after maintenance of flow-path components (e.g. pump seals, injection seal) to prove pressure tightness up to 600 bar.




Result Comment The testresult is ambiguos. Perform a Valve Test to verify that the pump has no leak.

Started at 15/01/2020 11:54:02 AM by LC-DAD\Tulip on LC-DAD using Agilent Lab Advisor - Version B.02.09 [195] - Basic

Results

Item	Value	Result
Channel	B	Done
Leakrate	µl/min	Failed
Current flow	0.000 ml/min	Done

Finished at 15/01/2020 11:59:55 AM

Signals**Limits**

Name	Limit
Maxiumum leak rate	3 µl/min

Signature: